


PHS MEDICAL<sup>®</sup>

C-port<sup>™</sup>



*Innovative solutions  
for long term vascular access...  
... from specialists!*

Single lumen C-Port™ Systems are available in the following sizes and configurations:


## Basic Cut-Down System

for implanting via venous cut-down\*\*

### ◆ Standard, venous

1 each - Port, Catheter, “Huber” Needle (straight), Blunt Needle, Vein Pick

Cath.-Size	Type	Ref.
5F	Polyurethane	SP-005CP
7.5F	Polyurethane	SP-075CP
8F	Silastic®-Silicone	SP-008CS
9.6F	Silastic®-Silicone	SP-096CS




\*appropriate sized introducer sets are available single, sterile for percutaneous implantation.

### ◆ Percutaneous Introducer Set

1 each: Split-Sheath Cannula w/Dilator, Introducer Needle, Guidewire with Straightener, 10ml Syringe

Size	for use with	Ref.
6F	Port w/ 5F Catheter	PI-0006
8F	Port w/ 7.5F Catheter	PI-0008
9F	Port w/ 8F Catheter	PI-0009
10F	Port w/ 9.6F Catheter	PI-0010



## Specifications, Port


Material:

Port: Titanium, MRI Compatible

Septum: Silicone elastomer


Weight: 12.5g

Volume: ~0.5ml



## Specifications, Catheter

Size	Material	O.D. (mm)	I.D. (mm)	Length (cm)
⊙ F5	Polyurethane	1.68	1.02	50
⊙ F7,5	Polyurethane	2.50	1.15	50
⊙ F8	Silicone	2.67	1.40	50
⊙ F9,6	Silicone	3.18	1.58	50



All catheters with tapered tip and length markings from 1 cm to 30 cm. Radiopaque.

**PHS MEDICAL®**  
**PHS Medical GmbH**

Leipzigerstr. 246, DE-34260 Kaufungen, Germany  
 Tel: +49(0)5605/94 96-0 Fax: +49(0)5605/94 96-99  
 e-mail: info@phs-medical.de