

CONDOR GOLD LINE Urology Surgical Retractor Systems

Precision, perfectly established processes and exactly precise work – this is what is required and necessary for urological interventions in operating theatres around the world. To ensure and to perfectly support this is a characteristic of the efficiency of CONDOR GOLD LINE Surgical Retractor Systems.

Today, these systems belong to the standard surgery equipment used in clinics throughout the world. And there is a good reason for it: **CONDOR Surgical Retractor Systems provide surgeons with the highest possible flexibility and maximum precision adjustment at the same time.**

This is offered on the basis of an intelligent and modular design system which in combination with retractor models used in special interventions can be exactly adjusted to the respective surgical procedure. Our surgical retractor systems can be combined individually and are fully compatible due to the modular design. You can choose freely which retractors you order from our different sets.

The following surgical retractor sets are available for the field of urology. Compose your personal surgery set.

Abdominal Maxi



Perineal Mini



Paediatric urology



All urology sets have one thing in common: they offer the core elements and advantages of the CONDOR Basic Set which allows for excellent settings and freedom of movement in all urological interventions by means of its unique ball elements. Each of the modular systems can be put to multifunctional use and can be perfectly positioned with just a few quick swift actions: patented blade holders ensure an optimum in-depth adjustment – tip-focused. One more advantage: through the integrated isolation in the central holder the risk for the patient is excluded by means of high-frequency leakage current.

You are convinced of the advantages? You would like to familiarise yourself with the Surgical Retractor System? We would be pleased to put a courtesy test equipment at your disposal.



Abdominal: Precision for abdominal urological interventions

Highest quality, easy handling and modular flexibility are the performance characteristics of the Abdominal Set. One more advantage: it can be put to multi-functional use. The set can be used in abdominal interventions such as, for example, cystectomy and radical prostatectomy as well as in nephrectomy when the patient is in a lateral position.

Flexible and three-dimensional site adjustments are thus possible in all urological interventions. And that is not all: special urological retractors form the basis for exact site adjustments in all surgical interventions. The bladder retractor comes in different sizes and allows for the tensing of the bladder in prostate surgeries.



An overview of the advantages of the Abdominal Set:


- Diverse utilisation, e.g. in cystectomy, radical prostatectomy and nephrectomy**
- Flexible and three-dimensional site adjustment**
- Flexible use of small and large retractors**
- Bladder retractor for tensing of the bladder**
- Tip-focused adjustment of retractors by means of patented blade holders**
- Individual lifting of the abdominal wall / costal arches**

Mini Urology Set

Mini is only the name of the set – the advantages provided to you by the set are absolutely maxi: the mini set convinces with its flexibility and perfect settings. Especially in perineal prostatectomy the mini system comes into operation. It is also designed for urological interventions in children.

Perineal Mini

Small retractors some of which are also elastic provide open access even in difficult site settings. The system is especially a first-choice option in a difficult positioning of patient. The flexible application of cold light which can be combined with the mini set is one of many examples for the fulfilment of our expectation to set the highest standards with regard to the development and further development of our systems.



An overview of the advantages of the Perineal Mini Set

- Tip-focused setting options**
- Special flexible retractors in different sizes**
- To be combined with cold light**
- Site adjustment in difficult positioning possible**



Paediatric urology

Urological interventions, especially in children, require maximum precision. Small incisions are required with regard to pyeloplasty or antireflux plasty. The adequate prerequisites are offered by the CONDOR Paediatric Urology Set. In addition to small elastic retractors, such as the Harrington retractor, the distance holder is part of the set's components. The central holder can be lowered to the body level by means of this distance holder. To the surgeon this means: perfect freedom of movement at the operating table. Mounting the system on the table guarantees absolute freedom of arm movement.




The Paediatric Urology Set can be perfectly combined with cold light. The site is ideally illuminated in depth, even with the smallest incisions.

The flexibility of the Paediatric Urology Set is especially demonstrated in pyeloplasty in a lateral position. With regard to this procedure, the surgery team benefits especially from the individually variable settings on either side of the ball joints: the surgical retractor system can ideally be adjusted to the patient's position.



An overview of the advantages of Paediatric Urology Set:

- Absolute freedom of movement and freedom of arm movement by means of the CONDOR Distance Holder
- Utilisation of small, flexible and elastic retractors
- Individual setting of each side by ball joints
- Can be combined with cold light



Familiarise yourself with the CONDOR Urology Surgical Retractor System. Use our offer of courtesy test equipment. A qualified sales person will assist you during first use. You will be convinced of the advantages and the easy handling!

Internationally recognised and utilised worldwide: CONDOR GOLD LINE



An overview of the advantages of the CONDOR Surgical Retractor System:

- Interchangeable on a modular basis
- Compatibility of retractors
- Multi-functionality
- 3-D setting
- Simple operability
- In-depth site adjustment
- Electrical safety
- User training through surgical accompaniment
- Easy cleaning of the device
- Optimisation of operational costs
- Certified in accordance with EN ISO 13485:2003 CE 0044
- Relief of qualified personnel
- Made in Germany

Apart from additional urology sets our portfolio also includes additional systems for gynaecology, endoprosthesis and general surgery.

The individual modular design of the system which can be pieced together on the platform of the basic version as required and in a fully compatible way meets the high expectations of customised solutions for urological surgeries.



CONDOR[®] GmbH
MEDICALTECHNIK

CONDOR GmbH MEDICALTECHNIK | Franz-Kleine-Str. 13 | 33154 Salzkotten | Germany

Telephone: +49 (0) 52 58 - 99 16-0 | Fax: +49 (0) 52 58 - 99 16-16 | Email: info@condor-med.de | www.condor-med.com

CONDOR[®] is a registered trademark of the CONDOR GmbH in Germany and in other countries. All rights reserved, particularly for utility models, designs and patents. Subject to technical modifications and delivery possibilities. The designations reproduced in this document can be brands in which the utilisation through third-parties for own purposes can infringe upon the proprietor's rights. Status: October 2007.