

# EMERGENCY WARMER ESH 04

VERSATILE. FOR MOBILE USE.

The Emergency warmer **ESH 04** was specially designed for mobile applications. It can warm up to two 500 ccm infusion bottles to 37°C and retains this temperature at the place of use for up to 90 minutes. With an integrated pressure cuff the emergency warmer can also be used as a stand-alone solution. It does not require an additional holder for carrying, the solutions are dispensed directly from the warmer.



## EMERGENCY WARMER ESH 04

- Designed for mobile use
- Versatile, compact size, low weight
- Fits into emergency kit
- Maintains the temperature for up to 90 minutes at the site of use
- 24 h continuous operation possible
- Integrated pressure cuff
- Solution dispensed directly from warmer

### TECHNICAL DATA:

<b>Device:</b>	Emergency warmer
<b>Type designation:</b>	ESH04
<b>Voltage:</b>	12V CC
<b>Power consumption:</b>	max. 20 W
<b>Type:</b>	mains and rechargeable battery operated
<b>Degree of protection against ingress of liquids:</b>	IPX splash proof
<b>Degree of protection against electric shock:</b>	B
<b>Dimensions:</b>	160 x 105 x 300 mm
<b>Weight:</b>	1.3 kg
<b>Classification:</b>	IIa according to rule 9
<b>Type of operation:</b>	continuous operation
<b>Mains adaptor</b>	
<b>Type designation:</b>	24W P2EFSW3
<b>Input:</b>	110-240V AC, 50-60 Hz, 550 mA
<b>Output:</b>	12V CC, 2A max.

ESH 04

E. BIEGLER GmbH  
Allhangstrasse 18a  
A-3001 Mauerbach (Vienna)  
Austria

Tel: +43 1 979 21 05  
Fax: +43 1 979 21 05 16  
office@biegler.com  
www.biegler.com

**BIEGLER**  
MEDIZIN ELEKTRONIK

CE  
0123

